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## Références

1. Duncan, H.F., Galler, K.M., Tomson, P.L. et al. (2019) 'European Society of Endodontology position statement: Management of deep caries and the exposed pulp', *International Endodontic Journal*, 52(7), pp. 923–934. Available at: [https://britishendodonticsociety.org.uk/\\_userfiles/pages/files/duncan\\_et\\_al2019international\\_endodontic\\_journal.pdf](https://britishendodonticsociety.org.uk/_userfiles/pages/files/duncan_et_al2019international_endodontic_journal.pdf)
2. Islam, R., Islam, M.R.R., Tanaka, T. et al. (2023) 'Direct pulp capping procedures – Evidence and practice', *Japanese Dental Science Review*, 59, pp. 48–61. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC9985044/>
3. Novotna, B., Holik, P., Morozova, Y. et al. (2024) 'Evaluation of Cytotoxicity of the Dental Materials TheraCal LC, TheraCal PT, ApaCal ART and Biodentine Used in Vital Pulp Therapy: In Vitro Study', *Dentistry Journal*, 12(8), p. 249. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC11352889/>
4. Bakhtiar, H., Nekoofar, M.H., Aminishakib, P. et al. (2017) 'Human Pulp Responses to Partial Pulpotomy Treatment with TheraCal as Compared with Biodentine and ProRoot MTA: A Clinical Trial', *Journal of Endodontics*, 43(11), pp. 1786–1791. Available at: <https://pubmed.ncbi.nlm.nih.gov/28822566/>
5. Lee, H., Shin, Y., Kim, S.O. et al. (2015) 'Comparative Study of Pulpal Responses to Pulpotomy with ProRoot MTA, RetroMTA, and TheraCal in Dogs' Teeth', *Journal of Endodontics*, 41(8), pp. 1317–1324. Available at: <https://pubmed.ncbi.nlm.nih.gov/26015158/>
6. Schuster, L., Sielker, S., Kleinheinz, J. & Dammaschke, T. (2025) 'Effect of light-cured pulp capping materials on human dental pulp cells in vitro', *International Endodontic Journal*, 58(7), pp. 1060–1072. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC12160968/>
7. Dureja A, Acharya SR, Kini S, Mayya A, Hedge V. Biocompatibility and Performance of Dental Composite Restorations: A Narrative Review on Free Monomer Release, Concerns and Solutions. *Engineering Proceedings*. 2023; 59(1):160. Available at: <https://doi.org/10.3390/engproc2023059160>
8. Hamid A, Hume WR, The effect of dentine thickness on diffusion of resin monomers in vitro. *Journal of Oral Rehabilitation*, 1997; 24: 20-25. Available at: <https://doi.org/10.1111/j.1365-2842.1997.tb00255.x>
9. Putzeys E, Corneliu Duca R, Coppens L, Vanoirbeek J, Godderis L, Van Meerbeek B, Van Landuyt KL. In-vitro transdental diffusion of monomers from adhesives,

Journal of Dentistry, 2028; Volume 75: 91-97. ISSN 0300-5712. Available at:  
<https://doi.org/10.1016/j.jdent.2018.05.023>.

10. About, I. (2022) *Biodentine: Properties and Clinical Applications*. Switzerland: Springer.31-50 URL: <https://link.springer.com/book/10.1007/978-3-030-80932-4>
11. Poggio, C., Ceci, M., Dagna, A., Beltrami, R., Colombo, M., Chiesa, M. (2015) 'In vitro cytotoxicity evaluation of different pulp capping materials: a comparative study', *Arh Hig Rada Toksikol.* 66(3), pp. 181-8. Available at:  
<https://pubmed.ncbi.nlm.nih.gov/26444338/>
12. Sahin, N., Saygili, S. & Akcay, M. (2021) 'Clinical, radiographic, and histological evaluation of three different pulp-capping materials in indirect pulp treatment of primary teeth: a randomized clinical trial', *Clinical Oral Investigations*, 25(6), pp. 3945–3955. Available at: <https://pubmed.ncbi.nlm.nih.gov/33404764/>
13. Camilleri, J. (2011) 'Characterization and hydration kinetics of tricalcium silicate cement for use as a dental biomaterial', *Dental Materials*, 27(8), pp. 836–844. Available at:  
<https://www.sciencedirect.com/science/article/abs/pii/S0109564111001151>
14. Camilleri, J. (n.d.) *Biodentine™: The dentine in a capsule or more*. Septodont Learning. Available at:  
<https://www.septodontlearning.co.uk/sites/default/files/Biodentine-The-dentine-in-a-capsule-or-more-by-Josette-Camilleri-compressed.pdf>
15. Elbanna, A., Atta, D., Sherief, D.I. (2022) 'In vitro bioactivity of newly introduced dual-cured resin-modified calcium silicate cement', *Dental Research Journal (Isfahan)*, 28(19), p. 1. Available at: <https://pubmed.ncbi.nlm.nih.gov/35308449/>
16. Gong, V., Franca, R. (2017) 'Nanoscale chemical surface characterization of four different types of pulp-capping materials', *Journal of Dentistry*, 58, pp. 11-18. Available at: <https://pubmed.ncbi.nlm.nih.gov/28043847/>
17. Nowicka, A., Wilk, G., Lipski, M., Kotecki, J., Buczkowska-Radlińska, J. (2015) 'Tomographic Evaluation of Reparative Dentin Formation after Direct Pulp Capping with Ca(OH)<sub>2</sub>, MTA, Biodentine, and Dentin Bonding System in Human Teeth', *J Endod.*, 41(8):1234-40. Available at: <https://pubmed.ncbi.nlm.nih.gov/26031301/>
18. Watson, T.F., Atmeh, A.R., Sajini, S. et al. (2014) 'Present and future of glass-ionomers and calcium-silicate cements as bioactive materials in dentistry: biophotonics-based interfacial analyses in health and disease', *Dental Materials*, 30(1), pp. 50–61. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC3885799/>

19. Kaur, M., Singh, H., Dhillon, J.S. et al. (2017) 'MTA versus Biodentine™: Review of literature with a comparative analysis', *Journal of Clinical and Diagnostic Research*, 11(8), pp. ZG01–ZG05. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC5620936/>
20. Sajini, S., Atmeh, A., Banerjee, A. et al. (2022) 'Glass-ionomer and calcium silicate-based cements interactions with human dentine in health and disease: Two-photon fluorescence microscopy and Raman spectroscopy analysis', *Dental Materials*, 38(11), pp. 1710–1720. Available at: <https://www.sciencedirect.com/science/article/abs/pii/S0109564122002536?via%3Dihub>
21. Atmeh, A.R., Chong, E.Z., Richard, G. et al. (2012) 'Dentin-cement interfacial interaction: calcium silicates and polyalkenoates', *Journal of Dentistry*, 91(5), pp. 454–459. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC4077527/>
22. Schneider LF and al. *Shrinkage Stresses Generated during Resin-Composite Applications: A Review*. J Dent Biomech. 2010. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC2951111/>
23. Gallo M, Abouelleil H, Chenal JM, Adrien J, et al. Polymerization shrinkage of resin-based composites for dental restorations: A digital volume correlation study. *Dental Materials*, 2019. Volume 35, Issue 11: 1654-1664. ISSN 0109-5641. Available at: <https://doi.org/10.1016/j.dental.2019.08.116>.
24. Kwon Y, Ferracane J, Lee IB. Effect of layering methods, composite type, and flowable liner on the polymerization shrinkage stress of light cured composites. *Dent Mater*. 2012;28(7):801-809. doi:10.1016/j.dental.2012.04.028. Available at: <https://pubmed.ncbi.nlm.nih.gov/22575738/>
25. Aksel H, Küçükkaya Eren S, Askerbeyli Örs S, Karaismailoğlu E. Surface and vertical dimensional changes of mineral trioxide aggregate and biodentine in different environmental conditions. *J Appl Oral Sci*. 2018 Dec 10;27:e20180093. doi: 10.1590/1678-7757-2018-0093. Available at: <https://pmc.ncbi.nlm.nih.gov/articles/PMC6296284/>
26. Abdelmegid, F. Y., Salama, F. S., Al-Mutairi, W. M., et al. (2017) 'Effect of different intermediary bases on microleakage of a nanocomposite interface in Class II box cavities in primary teeth: An in vitro study', *International Journal of Artificial Organs*, 40(2), pp. 82-87. Available at: <https://pubmed.ncbi.nlm.nih.gov/28222209/>
27. Septodont internal document (Biodentine Expert Board - September 11, 2023).